

# Maloney Vision Institute Intralase Centration Nomogram OD

Version 3

8/14/07

| To Center on the Pupil, if you<br>move the cursor: | Then Move:                       |
|--|----------------------------------|
| <b>11</b> or more clicks <i>Left</i>               | <i>recenter the suction ring</i> |
| <b>10</b> Clicks    Left                           | <b>11</b> Extra Clicks    Left   |
| <b>8</b> Clicks    Left                            | <b>10</b> Extra Clicks    Left   |
| <b>6</b> Clicks    Left                            | <b>9</b> Extra Clicks    Left    |
| <b>4</b> Clicks    Left                            | <b>8</b> Extra Clicks    Left    |
| <b>2</b> Clicks    Left                            | <b>7</b> Extra Clicks    Left    |
| <b>0</b> Clicks                                    | <b>6</b> Extra Clicks    Left    |
| <b>2</b> Clicks    Right                           | <b>5</b> Extra Clicks    Left    |
| <b>4</b> Clicks    Right                           | <b>4</b> Extra Clicks    Left    |
| <b>6</b> Clicks    Right                           | <b>3</b> Extra Clicks    Left    |
| <b>8</b> Clicks    Right                           | <b>2</b> Extra Clicks    Left    |
| <b>10</b> Clicks    Right                          | <b>1</b> Extra Click    Left     |
| <b>12</b> Clicks    Right                          | <b>0</b> Extra Clicks            |
| <b>14</b> Clicks    Right                          | <b>1</b> Extra Clicks    Right   |
| <b>16</b> Clicks    Right                          | <b>2</b> Extra Clicks    Right   |
| <b>17</b> Clicks    Right                          | <b>3</b> Extra Clicks    Right   |
| <b>18</b> Clicks    Right                          | <b>4</b> Extra Clicks    Right   |
| <b>19</b> or more clicks <i>Right</i>              | <i>recenter the suction ring</i> |

**Instructions:**

- 1) Count horizontal and diagonal clicks. Ignore vertical clicks.
- 2) Round down in-between values (e.g. 3 clicks left reads the row for 2 clicks left)
- 3) Copyright Robert K. Maloney, 2007. Permission is hereby given to any surgeon to reproduce this for personal use